


Alkyne-PEG4-NHS ester

<http://www.lumiprobe.com/p/alkyne-peg4-nhs-ester-1393330-40-9>

Alkyne-PEG4-NHS ester is a bifunctional molecule containing NHS ester and alkyne functional groups. NHS ester group is used for the labeling of amine groups, like lysine and terminal NH₂, in proteins and peptides. The terminal alkyne group is useful for copper-catalyzed click chemistry (CuAAC) or Sonogashira coupling.

PEG4 is a hydrophilic linker providing a good separation between the parts of the conjugate.



Structure of Alkyne-PEG4-NHS ester

General properties

Appearance: light yellow oil

Molecular weight: 401.41

CAS number: 1393330-40-9

Molecular formula: C₁₈H₂₇NO₉

IUPAC name: 4,7,10,13,16-Pentaoxanonadec-18-ynoic acid, 2,5-dioxo-1-pyrrolidinyl ester

Solubility: good in water and polar organic solvents

Quality control: NMR ¹H and HPLC-MS (95+%)

Storage conditions: 12 months after receipt at -20°C in the dark. Transportation: at room temperature for up to 3 weeks. Desiccate.

Legal statement: This Product is offered and sold for research purposes only. It has not been tested for safety and efficacy in food, drug, medical device, cosmetic, commercial or any other use. Supply does not express or imply authorization to use for any other purpose, including, without limitation, in vitro diagnostic purposes, in the manufacture of food or pharmaceutical products, in medical devices or in cosmetic products.