


Biotin-PEG4-azide

<http://www.lumiprobe.com/p/biotin-peg4-azide>

Biotin-PEG4-azide is a click chemistry biotinylation tool. This reagent allows the labeling of alkynylated molecules (such as DNA, oligonucleotides, and proteins) with biotin via copper-catalyzed or copper-free click reaction. Biotin-labeled biomolecules can be bound to avidin or streptavidin for further purification and detection.

This structure of biotin azide features a long hydrophilic PEG4 linker that separates biotin residue from the target molecule for efficient binding with streptavidin. Its linker also enhances aqueous solubility to facilitate conjugation.

The azide can be conjugated with various molecules using our recommended protocol.



Structure of Biotin-PEG4-azide

General properties

Appearance: beige solid

Molecular weight: 444.56

CAS number: 875770-34-6

Molecular formula: $C_{18}H_{32}N_6O_5S$

Solubility: in DMF, DMSO

Quality control: NMR 1H and HPLC-MS (95+%)

Storage conditions: 24 months after receipt at $-20^{\circ}C$ in the dark. Transportation: at room temperature for up to 3 weeks. Desiccate.