



JOE alkyne, 6-isomer

<http://www.lumiprobe.com/p/joe-alkyne-6>

JOE dye is a fluorophore similar to TET but with less hydrophobic properties. The dye is widely used for the labeling of nucleic acids. JOE alkyne can be conjugated with azides using a copper-catalyzed click chemistry reaction.



Structure of 6-JOE alkyne



Absorption and emission spectra of JOE dye

General properties

Appearance:	orange solid
Mass spec M+ increment:	541.0
Molecular weight:	542.32
Molecular formula:	C ₂₆ H ₁₇ NCl ₂ O ₈
Solubility:	good in DMF, DMSO
Quality control:	NMR ¹ H, HPLC-MS (95%)
Storage conditions:	Storage: 24 months after receipt at -20°C in the dark. Transportation: at room temperature for up to 3 weeks. Avoid prolonged exposure to light. Desiccate.
Legal statement:	This Product is offered and sold for research purposes only. It has not been tested for safety and efficacy in food, drug, medical device, cosmetic, commercial or any other use. Supply does not express or imply authorization to use for any other purpose, including, without limitation, in vitro diagnostic purposes, in the manufacture of food or pharmaceutical products, in medical devices or in cosmetic products.

Spectral properties

Excitation/absorption maximum, nm:	533
ε, L·mol ⁻¹ ·cm ⁻¹ :	75000
Emission maximum, nm:	554
Fluorescence quantum yield:	0.61
CF ₂₆₀ :	0.36
CF ₂₈₀ :	0.28