
QuDye dsDNA BR Reagent, 200×

<http://www.lumiprobe.com/p/qudye-dsdna-br-reagent>

QuDye dsDNA BR Reagent is highly selective to dsDNA over ssDNA, RNA, free nucleotides, or protein contaminants. The 200× formulation is also a component of the nucleic acid quantification [QuDye dsDNA BR Assay Kit](#) for dsDNA samples measurements using all popular readers.

General properties

Appearance: red-orange liquid

Quality control: NMR ¹H, HPLC-MS (95%)

Storage conditions: Storage: 12 months after receipt at -20°C in the dark. Transportation: at room temperature for up to 3 weeks. Avoid prolonged exposure to light. Desiccate.

Legal statement: This Product is offered and sold for research purposes only. It has not been tested for safety and efficacy in food, drug, medical device, cosmetic, commercial or any other use. Supply does not express or imply authorization to use for any other purpose, including, without limitation, in vitro diagnostic purposes, in the manufacture of food or pharmaceutical products, in medical devices or in cosmetic products.