


TBTA ligand

<http://www.lumiprobe.com/p/tbta-ligand>

TBTA (Tris(benzyltriazolylmethyl)amine) is a tertiary amine with three 1,2,3-triazole groups. It is one of the most widely used water-insoluble ligands for copper-catalyzed azide-alkyne cycloadditions (CuAAC). TBTA complexes with copper(I) and promotes click reaction catalysis through the stabilization of the copper(I)-oxidation state.



Structure of TBTA ligand

General properties

Appearance: white solid

Molecular weight: 530.64

CAS number: 510758-28-8

Molecular formula: $C_{30}H_{30}N_{10}$

IUPAC name: tris[(1-benzyl-1H-1,2,3-triazol-4-yl)methyl]amine

Solubility: good in DMSO, DMF, DCM, THF, Chloroform

Quality control: NMR 1H and HPLC-MS (95+%)

Storage conditions: 24 months after receipt at $-20^{\circ}C$ in the dark. Transportation: at room temperature for up to 3 weeks. Desiccate.

Legal statement: Product is offered and sold for research purposes only. Product is not tested for safety and efficacy in food, drug, medical device, cosmetic, no express or implied authorization to use for any other purpose, including, without limitation, in vitro diagnostic purposes, for humans or animals or for commercial purposes.